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PATIENT LEAFLET

1. General information

It is worth remembering that, although the materials used to make class III implantable medical devices are highly biocompatible and with a high safety profile, the success of surgery is largely influenced by compliance with post-operative indications and adherence to the rehabilitation protocol.

2. Medical device identification

Medical device identification name: DAC®.

In order to adapt to different surgical requirements and to allow a uniform coating of different surfaces, the product is offered on its own or as a kit. The product is available in different formats, quantities, sizes and physical state (powder or lyophilised product).

DAC® is a class III, implantable, single-use, bioresorbable medical device.

DAC® is a product designed for the preparation and application of a bio-absorbable hydrogel (gel whose liquid phase is water) coating for use in orthopaedic, dental and trauma surgery. The product is indicated in orthopaedic surgery in general including spine surgery, maxillofacial surgery, thoracic surgery, trauma surgery, joint surgery and in dentistry.

The device is indicated for:

- a) Patients of either sex, with mature bone age;
- b) Patients undergoing uncemented prosthesis implantation / joint arthroplasty;
- c) Patients undergoing open or closed reduction of skeletal fracture with respective internal or external fixation;
- d) Patients requiring internal fixation and stabilisation with bone-consolidating devices;
- e) Adult patients undergoing dental surgical implants;
- f) Patients undergoing spine surgery, maxillofacial surgery, thoracic surgery involving the skeletal system;
- g) Patients undergoing reconstructive tendon and ligament surgery involving the use of bone anchoring devices.

The function of this hydrogel is to create a temporary physical barrier on the surface of implantable devices against bacterial adhesion; it is therefore a preventive measure against bacterial adhesion, colonisation and biofilm formation (which causes bacterial infections) on the surface of implants. DAC® is a product designed for the prevention and not the treatment of infections.

The use of the medical device is restricted to specialised medical personnel. The product must be used in a sterile environment and/or respecting the criteria of hygiene, cleansing, cleaning normally required by Good Clinical Practice at the point of care.

3. Materials

The product consists of hyaluronic acid and polylactic acid. The biocompatible substances used in the manufacture are widely known molecules that have long been used in the production of implantable class III medical devices.

The product contains no medicinal substances, material of animal origin or blood components.

NOVAGENIT S.r.l. declares that phthalates, latex and nickel and other potentially hazardous materials are not present in the product. However, it cannot exclude the potential presence of technically unavoidable trace impurities due to the use of auxiliaries and additives in the production process.

4. Sensitivity to materials

Patients are encouraged to report suspected or established sensitivities to the components listed in section 3. of this document to their treating physician so that he can recommend the most appropriate treatment.

5. Use of diagnostic equipment

As bioresorbable materials, it is understood that Novagenit® products are compatible with the majority of MRIs found in hospitals.

Class III implantable Novagenit® medical devices fall under the category: "MR- safe". Such a category identifies products that can be freely introduced into MRI environments. An 'MR- safe' device is in fact made of electrically non-conductive, non-metallic and non-magnetic materials.

However, it will be the patient's responsibility to inform the facility that will perform the possible MRI examination, which will check for possible contraindications or effects that such materials may have on the examination.

6. Complications and undesirable effects

For the DAC® medical device there are no known complications or increased risks related to its use. The possibility of complications/adverse events may be related to inappropriate surgical technique, including incorrect positioning of the device, failure to follow post-operative indications or failure to follow the rehabilitation protocol.

Any unusual changes in the implant as well as falls or accidents should be reported, even if the device or surgical site did not appear to be damaged at the time.

7. Contraindications

There are no known contraindications. The hydrogel is resorbed within a few days, so it does not interfere with osseointegration processes. Toxicity to tissues in contact with the hydrogel has not been described.

The product should not be used in cases of hypersensitivity to any component (hyaluronic acid and/or polylactic acid).

Before proceeding to use the medical device, the problem must be assessed on a case-by-case basis, taking into account possible alternative surgical procedures.

8. Use of the device in special subjects

The effectiveness of the devices has not been established in the case of

- paediatric patients;
- pregnancy;
- lactation;
- renal diseases;
- neoplasia;
- cardio-vascular diseases;
- alcoholism;
- incapacitated patients unable to understand and follow the rehabilitation programme.

NOTE:

The complete list of information identifying the individual device, including the name, serial number, lot number, UDI and model of the device is indicated on the implant carrier card given to the patient by the Hospital Centre.

Any updates to the information in this document will be made available to the patient via the website www.novagenit.com.

www.novagenit.com

NOVAGENIT S.r.l.

Viale Trento 115/117
38017 – Mezzolombardo (TN)-ITALY
Tel. +39 0461 1916500
Fax +39 0461 1916591
E-mail info@novagenit.com