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E-mail info@novagenit.com

1. General information

It is worth remembering that, although the materials used to make class III implantable medical devices are highly biocompatible and with a high safety profile, the success of surgery is largely influenced by compliance with post-operative indications and adherence to the rehabilitation protocol.

2. Medical device identification

NovaBioGraft is a collagen-based matrix, also referred to as a 'scaffold' for tissue regeneration. The NovaBioGraft membrane is ready-to-use and is available in three product families: monophasic, biphasic and collagen-only. NovaBioGraft is available in different formats. Easily shaped and adapted to the defect site, it can also be applied arthroscopically.

NovaBioGraft® is a class III, implantable, single-use, bioresorbable medical device.

The device is intended to serve as a filler for bone and/or osteochondral and/or chondral defects and as a scaffold for tissue, chondral and/or osteochondral regeneration.

NovaBioGraft, when applied in a defect, guides the undifferentiated mesenchymal cell component to cartilage neoformation and bone regeneration.

In particular, the collagen-only version is indicated for tissue and chondral regeneration, the biphasic version (collagen+ HA-βTCP granules) is indicated for the regeneration of osteocartilaginous defects, and the monophasic version (collagen+ HA-βTCP granules) is indicated for bone defects.

The membrane is suitable for use in orthopaedic surgery and traumatology, but also in maxillo-facial surgery.

The use of the medical device is restricted to specialised medical personnel. The product must be used in a sterile environment and/or respecting the criteria of hygiene, cleansing, cleaning normally required by Good Clinical Practice at the point of care.

3. Materials

NovaBioGraft is composed of type I equine collagen. The collagen is extracted from tendons of equines, which come from controlled and licensed farms according to current European regulations. The collagen, extracted from tendons of healthy animals, is free of any transmissible virus or infectious agent. Equine tendons are free of BSE (Bovine Spongiform Encephalopathy) or risk of viral contamination.

The product contains no medicinal substances or blood components.

In the monophasic and biphasic versions, type I equine collagen is combined with beta-tricalcium phosphate hydroxyapatite HAβTCP granules. Both materials (collagen and HAβTCP granules) are widely known and used in medical devices.

NOVAGENIT S.r.l. declares that phthalates, latex and nickel and other potentially hazardous materials are not present in the product. However, it cannot exclude the potential presence of technically unavoidable trace impurities due to the use of auxiliaries and additives in the production process.

4. Sensitivity to materials

Patients are encouraged to report suspected or established sensitivities to the components listed in section 3. of this document to their treating physician so that he can recommend the most appropriate treatment.

5. Use of diagnostic equipment

As bioresorbable materials, it is understood that Novagenit® products are compatible with the majority of MRIs found in hospitals.

Class III implantable Novagenit® medical devices fall under the category: "MR- safe". Such a category identifies products that can be freely introduced into MRI environments. An 'MR- safe' device is in fact made of electrically non-conductive, non-metallic and non-magnetic materials.

However, it will be the patient's responsibility to inform the facility that will perform the possible MRI examination, which will check for possible contraindications or effects that such materials may have on the examination.

6. Complications and undesirable effects

There are no known complications or increased risks related to the use of the device.

Intraoperative and postoperative complications associated with the type of surgery may occur, such as infection, allergic reaction to the implant, pain, migration of the device, oedema, protracted swelling of the treated site, mild hyperpyrexia, joint stiffness, delayed healing of the defect. The same reactions can occur, as demonstrated in the literature, with the use of autologous bone.

These complications may require revision surgery with removal of the device.

The possibility of complications/adverse events may also be related to inappropriate surgical technique including incorrect device placement, failure to follow post-operative indications or failure to follow the rehabilitation protocol.

7. Contraindications

- Do not use the device in case of: presence of infection at the site to be treated, osteomyelitis or systemic infection;
- age over 60 years;
- osteochondral lesions greater than 10 cm²; obesity (BMI > 30);
- autoimmune diseases;
- diabetes;
- fractures secondary to metastases;
- osteoporotic fractures;
- vascular disorders;
- immunosuppression;
- presence of metabolic or systemic disorders that may interfere with bone healing processes;
- patients on steroid therapy;
- coagulation disorders; hypercalcaemia;
- established individual hypersensitivity to collagen type I of equine origin and calcium phosphate salts.

8. Use of the device in special subjects

The effectiveness of the devices has not been established in the case of

- paediatric patients;
- pregnancy;
- lactation;
- renal diseases;
- neoplasia;
- cardio-vascular diseases;
- alcoholism;
- incapacitated patients unable to understand and follow the rehabilitation programme.

NOTE:

The complete list of information identifying the individual device, including the name, serial number, lot number, UDI and model of the device is indicated on the implant carrier card given to the patient by the Hospital Centre.

Any updates to the information in this document will be made available to the patient via the website www.novagenit.com.

www.novagenit.com

NOVAGENIT S.r.l.

Viale Trento 115/117
38017 – Mezzolombardo (TN)-ITALY
Tel. +39 0461 1916500
Fax +39 0461 1916591
E-mail info@novagenit.com