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E-mail info@novagenit.com

1. General information

It is worth remembering that, although the materials used to make class III implantable medical devices are highly biocompatible and with a high safety profile, the success of surgery is largely influenced by compliance with post-operative indications and adherence to the rehabilitation protocol.

2. Medical device identification

The resorbable endorthesis for flat foot is a class III medical device, implantable, bioresorbable, disposable, designed to correct the flat foot deformity generated by excessive movement of the talus in relation to the calcaneus. Specifically, the endorthesis, appropriately inserted into the sinus tarsus, the natural cavity between the talus and calcaneus, limits its movement until correct talocalcaneal alignment is re-established; the endorthesis reabsorbs over time without the need for a second operation to remove it.

This device has the sole function of assisting bone welding and is not designed as a replacement for normal bone structures.

The resorbable endorthesis for flat feet is available in 2 sizes (8 and 10), with the dimensions, expressed in mm, referring to the diameter of the external cylindrical body (dowel). The choice of the right implant size for a given patient must be determined by taking into consideration the patient's height, weight, functional requirements and anatomical structure.

Indications for the use of this internal fixation device include the following malformation corrections:

- Flexible flatfoot in children from 8 years of age onwards, preferably from 8 to 12 years of age, considering that beyond the age of 12 years, the due accessory times are performed more frequently.
- Flat foot in adults
- Flat foot with tarsal fusion
- Congenital vertical talus
- Posterior tibial tendon dysfunction.

The use of the medical device is restricted to specialized medical personnel. The product must be used in a sterile environment and/or respecting the criteria of hygiene, cleansing, cleanliness normally required by Good Clinical Practice at the point of care. Implantation must be carried out using specific instrumentation.

3. Materials

The resorbable endorthesis for flat feet is made entirely (100%) of Poly (L-Lactide). The biocompatible substance used for manufacture is a molecule that is widely known and has long been used for the production of implantable class III medical devices.

The product contains no medicinal substances, material of animal origin or blood components.

NOVAGENIT S.r.l. declares that phthalates, latex and nickel and other potentially hazardous materials are not present in the product. However, it cannot exclude the potential presence of technically unavoidable trace impurities due to the use of auxiliaries and additives in the production process.

4. Sensitivity to materials

Patients are encouraged to report suspected or established sensitivities to the components listed in section 3. of this document to their treating physician so that he can recommend the most appropriate treatment.

5. Use of diagnostic equipment

As bioresorbable materials, it is understood that Novagenit® products are compatible with the majority of MRIs found in hospitals.

Class III implantable Novagenit® medical devices fall under the category: "MR- safe". Such a category identifies products that can be freely introduced into MRI environments. An 'MR- safe' device is in fact made of electrically non-conductive, non-metallic and non-magnetic materials.

However, it will be the patient's responsibility to inform the facility that will perform the possible MRI examination, which will check for possible contraindications or effects that such materials may have on the examination.

6. Complications and undesirable effects

In many cases, undesirable results are more likely to be attributable to the clinical situation than to the device itself. Listed below are the most frequent undesirable effects involving the use of internal fixation devices:

- These devices can break when subjected to overload. Internal fixation devices are load-sharing devices designed to support bony surfaces in a position that facilitates their bonding.
- Conditions attributable to poor bone formation can cause loosening, cracking or breakage of the device or premature loss of rigid fixation of the bone.
- Improper positioning may result in insufficient correction of the malformation and/or cracking or fracture of the device.
- Increased fibrous tissue around the surgical site due to improperly placed implants.
- Early or late infection, deep or superficial.
- Deep vein thrombosis.
- Risk of subclinical nerve damage as a result of surgical trauma.

Known cases of material sensitivity reactions in patients following surgical implantation are quite rare, however, their role is pending further clinical investigation.

7. Contraindications

Conditions presenting a greater risk of failure include:

- Any active infection or suspected latent infection as well as obvious localised inflammation in or near the affected area.
- Compromised vascularisation that could inhibit adequate blood supply to the surgical site.
- Compromised bone stock due to disease, infection or previous implantation that is unable to provide adequate support for the devices.
- Material sensitivity, verified or presumed.
- Obesity. A patient who is severely overweight or obese risks overloading the implant with possible damage.
- Use of an implant that could interfere with anatomical structures or physiological performance.

- Any mental or neuromuscular disorder that could imply an unacceptable risk of loss of fixation or complications during postoperative care.
- Other medical or surgical conditions that might preclude the potential benefits of surgery.

8. Use of the device in special subjects

The effectiveness of the devices has not been established in the case of

- paediatric patients (age < 8 years);
- pregnancy;
- lactation;
- renal diseases;
- neoplasia;
- cardio-vascular diseases;
- alcoholism;
- incapacitated patients unable to understand and follow the rehabilitation programme.

NOTE:

The complete list of information identifying the individual device, including the name, lot number, UDI and model of the device is indicated on the implant carrier card given to the patient by the Hospital Centre.

Any updates to the information in this document will be made available to the patient via the website www.novagenit.com.

www.novagenit.com

NOVAGENIT S.r.l.

Viale Trento 115/117
38017 – Mezzolombardo (TN)-ITALY
Tel. +39 0461 1916500
Fax +39 0461 1916591
E-mail info@novagenit.com