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1. General information

It is worth remembering that, although the materials used to make class III implantable medical devices are highly biocompatible and with a high safety profile, the success of surgery is largely influenced by compliance with post-operative indications and adherence to the rehabilitation protocol.

2. Medical device identification

Medical device identification name: ElastiCo®.

ElastiCo is a Class III, implantable, single-use, bioresorbable medical device.

The ready-to-use porous ElastiCo membrane (commonly referred to as scaffold) is available in two configurations: three-dimensional sponge and two-dimensional film for use in the regenerative approach to the treatment of cartilage and tendon injuries. ElastiCo also acts as a barrier against peritendinous adhesions.

It serves as a structure for undifferentiated mesenchymal stromal cells from the blood and resorbs as the connective tissue neoforms and the defect heals completely.

The three-dimensional conformation of ElastiCo is indicated as a support for cartilage regeneration in cases of cartilage injuries of traumatic, post-traumatic, degenerative origin.

The time period of total replacement with new cartilaginous connective tissue depends on a number of variables such as quality and vitality of the recipient subchondral bone, volume and extent of the lesion, age of the patient, stabilisation of the defect, expected loading condition as well as other factors varying from patient to patient.

The two-dimensional shape of ElastiCo is indicated as a support for the regeneration of tendon connectives. The device in the two-dimensional form also acts as a barrier in the prevention of peritendinous fibroadherences.

In particular, it is indicated in:

- treatment of partial rotator cuff injuries of degenerative or traumatic origin that require surgery.
- treatment of tendon injuries of the Achilles tendon.
- treatment of patellar, biceps and triceps tendon injuries.
- prevention of peritendinous adhesions.

The use of the medical device is restricted to specialised medical personnel. The product must be used in a sterile environment and/or respecting the criteria of hygiene, cleansing, cleanliness normally required by Good Clinical Practice at the point of care.

3. Materials

The device is composed of type I equine collagen.

The collagen is extracted from tendons of equines, which come from controlled and licensed farms according to current European regulations.

The collagen, extracted from tendons of healthy animals, is free of any transmissible virus or infectious agent. Equine tendons are free of BSE (Bovine Spongiform Encephalopathy) or risk of viral contamination.

The product contains no medicinal substances or blood components.

NOVAGENIT S.r.l. declares that phthalates, latex and nickel and other potentially hazardous materials are not present in the product. However, it cannot exclude the potential presence of technically unavoidable trace impurities due to the use of auxiliaries and additives in the production process.

4. Sensitivity to materials

Patients are encouraged to report suspected or established sensitivities to the components listed in section 3. of this document to their treating physician so that he can recommend the most appropriate treatment.

5. Use of diagnostic equipment

As bioresorbable materials, it is understood that Novagenit® products are compatible with the majority of MRIs found in hospitals.

Class III implantable Novagenit® medical devices fall under the category: "MR- safe". Such a category identifies products that can be freely introduced into MRI environments. An 'MR- safe' device is in fact made of electrically non-conductive, non-metallic and non-magnetic materials.

However, it will be the patient's responsibility to inform the facility that will perform the possible MRI examination, which will check for possible contraindications or effects that such materials may have on the examination.

6. Complications and undesirable effects

For ElastiCo, there are no known complications or increased risks related to the use of the device.

Intraoperative and postoperative complications associated with the type of surgery may occur, such as infection, inflammation, allergic reaction to the implant, pain, migration of the device, oedema, protracted swelling of the treated site, mild hyperpyrexia, joint stiffness, delayed healing of the defect. The same reactions can occur, as demonstrated in the literature, with techniques commonly used for the treatment of cartilage defects or tendon injuries.

These complications may require revision surgery with removal of the device.

The possibility of complications/adverse events may also be related to inappropriate surgical technique including incorrect device placement, failure to follow post-operative indications or failure to follow the rehabilitation protocol.

7. Contraindications

Contraindications in the treatment of cartilage lesions:

- presence of infection at the site to be treated or systemic infection;
- age over 60 years;
- chondral lesions larger than 10 cm²; obesity (BMI > 30);
- autoimmune diseases;
- diabetes;
- vascular disorders;

- immunosuppression;
- metastases;
- presence of metabolic or systemic disorders that may interfere with cartilage healing processes;
- patients on steroid therapy;
- coagulation disorders;
- established individual hypersensitivity to collagen type I of equine origin.

Contraindications in tendon surgery:

- presence of active infection at the site to be treated or systemic infection;
- absence or poor vascularisation in the case of rotator cuff reconstruction;
- pathological soft tissue conditions that impair tendon fixation in the case of rotator cuff reconstruction;
- autoimmune pathologies;
- vascular disorders;
- presence of metabolic or systemic disorders that may interfere with tendon healing processes;
- coagulation disorders;
- established individual hypersensitivity to collagen type I of equine origin.

8. Use of the device in special subjects

The effectiveness of the devices has not been established in the case of

- paediatric patients;
- pregnancy;
- lactation;
- renal diseases;
- neoplasia;
- cardio-vascular diseases;
- alcoholism;
- incapacitated patients unable to understand and follow the rehabilitation programme.

NOTE:

The complete list of information identifying the individual device, including the name, serial number, lot number, UDI and model of the device is indicated on the implant carrier card given to the patient by the Hospital Centre.

Any updates to the information in this document will be made available to the patient via the website www.novagenit.com.

www.novagenit.com

NOVAGENIT S.r.l.

Viale Trento 115/117

38017 – Mezzolombardo (TN)-ITALY

Tel. +39 0461 1916500

Fax +39 0461 1916591

E-mail info@novagenit.com