

MATERIALS COMPOSING NOVAGENIT® IMPLANTS
MATERIALI CHE COMPONGONO GLI IMPIANTI NOVAGENIT®
MATÉRIAUX COMPOSANT LES IMPLANTS NOVAGENIT®
MATERIALIEN AUS NOVAGENIT® IMPLANTATEN
MATERIALES QUE COMPONEN LOS IMPLANTES NOVAGENIT®

REF	Note	Materials	Standard
NBS020600	Kit NovaBarrier S - 1x_60mg_1ml antiadhesive hydrogel (Extremity)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS021200	KIT NovaBarrier S-1x -120mg_2ml antiadhesive hydrogel (Extremity)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS033000	Kit NovaBarrier S- 300 mg_5ml antiadhesive hydrogel (vertebral)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS031800	Kit NovaBarrier S – 180 mg_3ml antiadhesive hydrogel (vertebral)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBL3210	NovaBarrier L- antiadhesive sheet 50x50x0,2mm-300mg (Extremity)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBL3300	NovaBarrier L – antiadhesive membrane 50x50x3mm-300mg (vertebral)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBL6300	NovaBarrier L – antiadhesive membrane 50x50x8mm-600mg (vertebral)	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS020300	Kit NovaBarrier S_1x_30mg_0,5ml_antiadhesive hydrogel_extremity	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS021220	Kit NovaBarrier S_1x_120 mg_2ml antiadhesive hydrogel_cannula_extremity	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS033020	Kit NovaBarrier S_1x_300 mg-5ml antiadhesive hydrogel_cannula_vertebral	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series

NBS031820	Kit NovaBarrier S_1x_180mg_3 ml antiadhesive hydrogel_cannula_vertebral	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS020620	Kit NovaBarrier S_1x_60mg_1ml antiadhesive hydrogel_cannula_extremity	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
NBS020320	Kit NovaBarrier S_1x_30mg_0,5ml natiadhesive hydrogel_cannula_extremity	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC010600	DAC - 1 ml. – dental	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC011200	DAC - 2 ml. – dental	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC001200	DAC - 2 ml. - ORTO	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC001800	DAC - 3 ml. - ORTO	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC003000	DAC - 5 ml. - ORTO	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC003002	DAC - 10 ml. - ORTO	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC003003	DAC - 15 ml. - ORTO	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC031800	DAC - 3 ml. - VERT	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC033000	DAC - 5 ml. - VERT	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC033002	DAC - 10 ml. - VERT	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC041800	DAC - 3 ml. - MAX	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC043000	DAC - 5 ml. - MAX	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC043002	DAC - 10 ml. - MAX	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series

DAC053000	DAC - 5 ml. - TOR	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC053002	DAC - 10 ml. - TOR	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DAC053003	DAC - 15 ml. - TOR	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DACL3000	DAC LYO_5 x 5 x 0,4 cm_ortho	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DACL6000	DAC LYO_5 x 5 x 0,8 cm_ortho	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DACL9000	DAC LYO_5 x 5 x 1,2 cm_ortho	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
DACL3010	DAC LYO_5 x 5 x 0,1 cm_ortho	hyaluronic acid+ poly D-L lactide	ASTM F 1635, ISO 13781, ASTM F 1925, ASTM F 2579, ISO 10993 series
ELP404010	ElastiCo sheet 40x40x0,1mm (Tendon)	type I equin collagen	ASTM F2212, ISO 22442-1/2/3
ELM404030	Elastico membrane/ElastiCo Sponge 40x40x3mm (Chondral)	type I equin collagen	ASTM F2212, ISO 22442-1/2/3
ELP404100	ElastiCo Film 40x40x1mm (Tendon)	type I equin collagen	ASTM F2212, ISO 22442-1/2/3
COL404030	NovaBiograft C 40x40x3mm (chondral)	type I equin collagen	ASTM F2212, ISO 22442-1/2/3
CHB404040	NovaBiograft OC 40x40x4mm (Osteocondral)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03
CHM404070	NovaBiograft B 40x40x7mm (Bone)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03
CHB404070	NovaBiograft OC 40x40x7mm (Osteocondral)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03

CHB404100	NovaBiograft OC 40x40x10mm (Osteocondral)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03
COL404070	NovaBiograft C 40x40x7mm (Chondral)	type I equin collagen	ASTM F2212, ISO 22442-1/2/3
CHM404040	NovaBiograft B 40x40x4mm (Bone)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03
CHM303040	NovaBiograft B 30x30x4mm (Bone)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03
CHM202040	NovaBiograft B 20x20x4mm (Bone)	type I equin collagen; HA-bTCP	ASTM F2212, ISO 22442-1/2/3, ISO 13779-3, ASTM F1088 – 18, ASTM F1185-03
2100008	BFFI tg 8	Poly L lactide	ISO 13781, ASTM F 1925, ASTM F 2579
2100010	BFFI tg 10	Poly L lactide	ISO 13781, ASTM F 1925, ASTM F 2579